


CardioNet

A new technology
eliminating distances
shortening the time
standardizing quality





The **CardioNet** software from **Cardios** is a new concept for transmission and reception of Holter recordings and reports with the **CardioLight** through the Internet or Intranet.

It's composed of two modules: **CardioNet Client** and **CardioNet Server**.

The **CardioNet Client** is installed at places where the hook-up is done (Generating Local Units), that may be clinics, doctor's offices or hospitals and allows, among other things, to transmit a recording through the internet or intranet to a **Central Analysis Unit**.

The **CardioNet Server** is installed at the Central Analysis Unit which may be located at hospitals, clinics, diagnostic centers or laboratories. It is integrated with the **CardioSmart** software and among other things, allows to:

- Perform all **CardioNet Client** functions.
- Receive **Holter** recordings from various **Local Units**.
- Transmit the reports to the **Local Units** that have sent the recording.

This makes possible to fully integrate the Local Units with the Central Analysis Unit regardless of the location, even at different cities, with great advantages for doctors and patients.

Important features:


1. Discards the traffic of tapes, memory cards or reports between the Local Units and the Central Analysis Units, increasing the security of the patient's data.



2. Eliminates costs of transport (mail, carrier, etc) that would be passed to the patient or would be assumed by the doctor.

3. Avoids the necessity of investment in an Analysis Central for doctors with small volume of recordings each month.

This new technology allows for any patient to enjoy the benefits of the same quality standard of his Holter exam done by Cardios' equipments, disregarding where he is located, with the same speed, trust and security that has when the tests are made at any Central Analysis Unit chosen by the doctor.



Cardio Sistemas Indl. Coml. Ltda
 Av. Paulista, 509 1st floor 01311-910 São Paulo Brasil
 Tel. 55 11 3883-3000 Fax 55 11 3883-3060
 www.cardios.com.br

